floorline LTC

ABCOGROUP OFFICE Solutions. Simplified.









Our Floorline LTC is the optimal bed solution for home care use or long-term care facilities. Its one piece design makes it easy to fold, store and transport. The Floorline LTC comes in two convenient options including one with a wood-grain finish that compliments any room.

Clinical and Patient Benefits

The Floorline LTC lowers to a distance of 9.5cm/3.7" from the floor. This distance, although seemingly small, has enormous benefits to both the patient and the facility. By just adding 5cm to the distance between the bed and the floor, the impact force of a fall is increased by 50%. This can be the difference between a simple bruise or a fractured hip for the patient that may result in an extented hospital stay or even death. In addition to the physical effects, fall injuries also have psychological effects on the patient as they experience a loss of confidence and fear after a fall.

Optimal Features

The Floorline LTC features standard positions

normally found only on hospital beds, such as Trendelenburg, reverse Trendelenburg and cardiac positioning. With a full assortment of accessories, an easy-to-use handset, and durable LINAK electronics, the Floorline LTC is the best choice for patient care.

Quality

Our products are designed and manufactured to meet all the requirements of relevant standards. At every stage of the development process, every single product is given extensive attention in design, development and manufacturing. This means that a user of this product feels confident with its reliability and the caregiver feels confident with its safety.

Accessories

ന

IV Pole

Adjustable height. Can be mounted on either end of bed. Maximum SWL 7kg/ 15.4lb. Part no. IVP



Care Assist Rail Designed to assist patients safely in and

out of a bed. Ergonomic design. Part no. 90702



Extension Kit Extends bed by 17.5cm/7" for taller patients. Part no. 90709



Handset Five function - included with bed. Part no. 75746 Handset holder Part no. LC2250, Lock Out Key Part no. 75790

90600 (EU), 90602 (UK),

90604 (NA), 90608 (AU)

9.5cm - 70cm / 3.7" - 27.5"

215cm x 90.5cm /

84

Technical Specifications

Floorline LTC

Part no. White Finish

Platform Position

Platform Dimensions (LxW)

Sleeping Dimensions (LxW)

(Min-Max)

Output Voltage

Input Voltage

Duty Cycle

Platform Panel Angles







Oxygen Holder

Self Help Pole

SHP600

Assists patient moving

within confines of

bed. Maximum SWL

75kg/165.3lb. Part no.

Carries C size tank. Attaches to legs. Part no. 90707

90601 (EU), 90603 (UK),

90605 (NA), 90609 (AU)

215cm - 232.5cm /

232.5cm x 90.5cm /

84.5" - 91.5"

Features

Compact

Folds compactly for storage or transport.



Standard Positions Includes auto-contour, backrest, kneebreak, Trendelenburg, reverse Trendelenburg, cardiac.

Premium Casters

One central locking and three steering/directional casters for premium maneuvrability.



Manufactured to all relevant standards including IEC 60601-2-52:2009, EN1970:2002, AS/NZS 3200.2.38:2007

84.5" x 35.5"	Extension Kit (L x W)	91.5" x 35.5"
200cm x 90cm / 78.7" x 35.5"	Sleeping Dimensions with Extension Kit (LxW)	217.5cm x 90cm / 85.5" x 35.5"
Max. 24V DC	IP Rating	IPX4
120V, 60Hz (90404 to 90607), else 230V, 50Hz	Product Weight	106kg/234lb
Backrest 70°, Thigh 45°, Knee 110°, Calf 25°, Trendelenburg/ Reverse Trendelenburg 18°	Safe Working Load (SWL)	200kg/440lb
2 minutes operation, 18 minutes rest	Maximum User Weight	165kg/363lb
This product is propulational		

Part no. Woodgrain Finish

Overall Bed Length (Standard to

Extended with Extension Kit)

Platform Dimensions with

This product is manufactured in accordance with all relevant ISO standards.



Learn More About This Product



abco@abcogroup.ca